

Mepilex® Border Flex | Time to change

When to change dressing according to saturation

The ideal wound dressing is designed to minimize the level of tissue disturbance. Too frequent dressing changes can lead to reduced temperature of the wound bed and disturbance to healing cells.¹ The added benefit to an institution is a reduction in spend due to dressings being changed too often.

Saturation Levels:




0

Strikethrough at 0 edges.
Can keep in place.



1

Strikethrough at 1 edge.
Can keep in place.



2

Strikethrough at 2 edges.
Can keep in place.



3

The dressing is saturated when the strikethrough has reached 3 edges.




Mepilex® Border Flex offers up to a 7 day wear time.


Mepilex Border Flex is designed for superior^{2,3} fluid retention. **Initial strikethrough does not necessarily mean that the dressing is saturated.** The dressing was constructed to make the strikethrough visible so that you can determine when it is time to change it.

Mepilex Border Flex may be left in place for 7 days depending on the condition of wound and surrounding skin.

Designed for undisturbed wound healing



Proven to handle more fluid than other dressings:



Fluid handling capacity is the amount of fluid absorbed by the dressing and its moisture loss through the backing film. A recent laboratory test showed that Mepilex Border Flex had superior fluid handling capacity to other brands.⁶

References: 1. Taking the trauma out of wound care: the importance of undisturbed healing by M. Rippon et al Journal of Wound Care Vol 21, No 8, August 2012. 2. White R. A Multinational survey of the assessment of pain when removing dressings. Wounds UK, 2008. 3. White R. et al. Evidence for atraumatic soft silicone wound dressing use. Wounds UK, 2005. 4. Mepilex Border Flex –Waterproofness. Report no. PD-532095_01. Data on file. 5. Viral Penetration Resistance. Report no. PD-535090_02. Data on file. 6. Mepilex Border Flex – Fluid handling capacity. Report no. PD-527642_01. Data on file. 7. Mepilex Border Flex—Bacteria encapsulation. Report no. PD-537072_04. Data on file. 8. Mepilex Border Flex—Absorption and dispersion tests on inclined plane, viscous test solution. Report no. PD-528871_02. Data on file. 9. Wiberg A.B. et al. Preventing maceration with a soft silicone dressing: in-vitro evaluations. Poster presented at the 3rd Congress of the WUWHs, Toronto, Canada, 2008. 10. RISE Report, Fluid Handling Capacity, 2019. 11. Molnlycke Health Care. Allevyn Life Testing. Lab Report No. 20171024-007.

5550 Peachtree Parkway, Suite 500, Norcross, GA 30092 | 1-800-843-8497 | www.molnlycke.us/MBFlex

The Molnlycke, Mepilex and Safetac trademarks, names and logo types are registered globally to one or more of the Molnlycke Health Care Group of Companies. Distributed by Molnlycke Health Care US, LLC, Norcross, Georgia 30092. © 2020 Molnlycke Health Care AB. All rights reserved. 1-800-843-8497. MHC-2018-T00018 Aquacel is a registered trademark of ConvaTec. | Allevyn[®] is a registered trademark of Smith & Nephew, Inc. | Optifoam[®] is a registered trademark of of Medline Industries, Inc.