

# Bruker's Workflow and Recommendations for Flow Sort Stains


The following workflow has been used in generating single-cell functional data and serves as a guideline for flow sort staining with the Bruker system.

## In this Technical Note we outline:

- Workflow recommendations for stain combinations
- Bruker's standard protocol and flow sorting protocol



## Bruker's Standard Protocol and Recommended Flow Sorting Protocol



**Figure 1 | a) Standard Protocol.** Sample enrichment, cell stimulation, cell staining, load cells.  
**b) Flow Sorting Protocol.** Flow sorting, cell stimulation, cell staining, load cells.

## Prep, Run, Analyze

### Workflow Recommendations for Stain Combinations

Bruker has validated flow sort staining for CD4<sup>+</sup> and CD8<sup>+</sup> T cells on the IsoLight. Users can now flow sort cells prior to loading into their IsoLight system (Figure 1). The following workflow is recommended:

- 1. Thaw PBMCs.** Thaw cryopreserved PBMCs and resuspend cells with complete RPMI and IL-2.
- 2. Recover PBMCs.** Incubate cells overnight.
- 3. Stain PBMCs.** Using the recommended flow stains, as shown in Figure 1, stain the cell subset of interest CD4<sup>+</sup> and/or CD8<sup>+</sup> T cells. Use flow stains per manufacturer's directions.
- 4. Flow sort.** Sort the cells using standard flow sorting protocol.
- 5. Stimulate cells.** Collect sorted cell subset of CD4<sup>+</sup> and/or CD8<sup>+</sup> T cells. Stimulate the cells on plate for 1.5 hours with Cell Activation Cocktail (PMA and Ionomycin).
- 6. Thaw IsoCode chips.** Allow IsoCode chips to thaw for 30-60 minutes prior to use. Set up all liquid reagents in instrument.
- 7. Stain CD4<sup>+</sup>/CD8<sup>+</sup> T cells.** Using Bruker's cell stain 405, stain each cell subset. It is critical to use an Bruker validated stain after flow sorting.
- 8. Load cells onto chips.** Resuspend cells using complete RPMI to a density of  $1 \times 10^6$  cells/mL and load 30  $\mu$ L of cell suspension onto IsoCode chip.
- 9. Load all chips onto instrument and run experiment.**
- 10. Analyze data using IsoSpeak.** In IsoSpeak, label flow sort stains as red and cell stain 405 as violet.