



NovaCardio



Fully integrated cardiology solution

NovaCardio - part of the NovaPACS/Ncompass Enterprise Imaging System - marries a best-in-class multi-modality motion cardiac imaging viewer with quick, easy, accurate and accreditation-compliant reporting.

Cardiologist designed, web-based and scalable, your team can access NovaCardio from a dedicated workstation or any computer within the facility or remotely.

Novarad's user interface includes familiar icons and navigation tools to streamline the workflow. Now, physicians can maximize their user experience with intuitive, efficient interaction.

NovaCardio enables facilities to:

- Get up and running quickly and easily – our US-based physicians, clinical experts, technical staff and support teams will work to configure and maintain the system for you. Even more rapid deployment is possible for existing NovaPACS sites
- Streamline and customize workflow, speed information access, and improve report turnaround time (TAT)
- Eliminate manual /redundant data entry
 - Demographic /historical/blood lab data can flow from the EMR /modality directly to the report
 - Measurements transfer automatically from modality to report and drive report creation
- Increase profitability and lower costs through improved billing cycle time and compliance
- Stay compliant with standards organizations like IAC while producing concise yet comprehensive reports
- View/review both radiology and cardiology images in a single program
- Deploy our solution on-site or in the cloud




Cardiac study types include:

- Echocardiography (Adult, Pediatric, Fetal, Stress, TEE)
- Nuclear Stress Testing
- Vascular Ultrasound
- Electrophysiology
- Cardiac Catheterization
- Cardiac CT & MRI
- ECG




Features

- Diagnostic quality DICOM Image review allows side-by-side review of images and reports. Still frame images from studies can be transferred into final reports
- Modality-specific rules-based error checking and billing compliance ensures the highest quality reports and when data is missing or conflicts
- Smart technology assists in selecting the correct ICD-10 and CPT codes
- Reports feature drawings and automatically generated graphics
- Completed reports are immediately available to referring physicians via auto-fax and can be sent to your existing EMR/RIS/HIS
- Web-based and vendor-neutral reference viewer accessible in the EMR
- Smart worklist features can enforce contractual agreements (SLAs) and clinical protocols
- Secure, with both redundancy and offsite disaster recovery included
- Dashboards and Reports include lab and staff activity information and turnaround time
- Access archived studies instantly



Learn more at

www.novarad.net/novacardio



NOVARAD CORPORATE HEADQUARTERS

3152 N University Ave, Suite 250
Provo, Utah 84604
(877) 668-2723

UNITED KINGDOM

12 Kingsbury Trading Estate Church Ln.
Kingsbury London United Kingdom
NW9 8AU
Phone: +44 (0) 208 205 9500
Fax: +44 (0) 208 205 0585

LATIN AMERICA

2 Calle A 6-28 zona 10 Edificio
Verona, Oficina 502
Edificio Verona, Oficina 502
Guatemala, Guatemala

ASIA PACIFIC

4/F SycipLaw Center
105 Paseo De Roxas
San Lorenzo Village
Makati City, Philippines