







## MICROVAP Parts and Accessories List



*\*\*Please contact us for part numbers for acid resistant units\*\**

Item Description	Part Number	Picture
<p>Stainless Steel Needles, 4 in x 19 gauge, 304 grade shaft, brass plated hub</p> <p>Sold in sets of 12 or 50</p>	<p>Set of 12 NA0603</p> <p>Set of 50 NA1303</p>	
<p>FEP coated needles, 4 in x 19 gauge, for corrosive solvents</p> <p>Sold in sets of 12 or 50</p>	<p>Set of 12 NA0603-T</p> <p>Set of 50 NA1303-T</p>	
<p>Pasteur pipette adapter and flow controller, male luer to 1/8 hose barb</p> <p>Sold in sets of 12</p>	<p>Set of 12 NA0637</p>	
<p>Plastic luer plug for unused luer fittings</p> <p>Sold individually</p>	<p>NA1307</p>	
<p>Connector tube for gas supply, 6 ft coiled polyurethane plastic</p> <p>Sold individually</p>	<p>NA1101</p>	

# Organomation®

Item Description	Part Number	Picture
<p>Luer and washer 10/32 in, chrome plated brass for MICROVAP nitrogen distribution manifold</p> <p>Sold individually</p>	P0607	
<p>Pressure reducing regulator for gas, 30-150 psig in, 0-30 psig out</p>	NA0630	
<p>Nitrogen Gas Generator</p> <p>20 LPM up to 48 sample positions</p>	NA1955	
<p>Custom sized inserts for 15 Position MICROVAP</p> <p>Sold individually or in sets of 15</p>	<p>Individuals NA1801</p> <p>Set of 15 NA1802</p>	
<p>Customized 24 position heat block</p> <p>Sold individually</p>	NA1824	
<p>15 Position Conversion Kit</p> <p><b>Includes:</b></p> <ul style="list-style-type: none"> <li>• 15 position heat block</li> <li>• set of 15 customized inserts</li> <li>• 15 position manifold</li> </ul>	<p>Heat Block: NA1815</p> <p>Inserts: NA1802</p> <p>Manifold: NA1805</p>	

# Organomation®

Item Description	Part Number	Picture
<p>24 Position Conversion Kit</p> <p><b>Includes:</b></p> <ul style="list-style-type: none"><li>• 24 position heat block</li><li>• 24 position manifold</li></ul>	<p>Heat Block: NA1824</p> <p>Manifold: NA1806</p>	 The image shows two components of the 24-position conversion kit. On the left is a rectangular, dark-colored heat block with a grid of 24 circular ports. On the right is a manifold assembly, which is a horizontal bar with 24 vertical tubes extending downwards, each ending in a small, dark-colored tip. The manifold is connected to a larger, more complex assembly at the top right.
<p>Thermometer adapter for bath calibration</p> <p>Compatible with single position and triple position MICROVAPs</p>	<p>NA1899</p>	 The image shows a single, conical, metallic thermometer adapter. It has a wide, flat base and tapers to a narrower top. The surface is polished and reflective.