












N-EVAP Parts and Accessories List







Please contact us for part numbers for acid resistant units

Item Description	Part Number	Picture
<p>Stainless Steel Needles, 4 in x 19 gauge, 304 grade shaft, brass plated hub</p> <p>Sold in sets of 12 or 50</p>	<p>Set of 12 NA0603</p> <p>Set of 50 NA1303</p>	
<p>Stainless Steel Needles, 4 in x 19 gauge, 304 grade shaft, plastic hub</p> <p>Sold in sets of 12 or 50</p>	<p>Set of 12 NA0604-19</p> <p>Set of 50 NA1304-19</p>	
<p>FEP coated needles, 4 in x 19 gauge, for corrosive solvents</p> <p>Sold in sets of 12 or 50</p>	<p>Set of 12 NA0603-T</p> <p>Set of 50 NA1303-T</p>	
<p>Nitrogen gas generator</p> <p>20 LPM up to 48 sample positions</p>	<p>NA1955</p>	
<p>Pasteur pipette adapter and flow controller, male luer to fluted nut</p> <p>Sold in sets of 12</p>	<p>Set of 12 NA0636</p>	







Organomation®

Item Description	Part Number	Picture
<p>Pasteur pipette adapter and flow controller, male luer to 1/8 hose barb</p> <p>Sold in sets of 12</p>	<p>Set of 12 NA0637</p>	
<p>Nylon 3-part tube fitting for N-EVAP valve tube</p> <p>Sold individually</p>	<p>P0634</p>	
<p>Luer and washer in chrome plated brass, for N-EVAP valve tube</p> <p>Sold individually</p>	<p>P0607</p>	
<p>Valve tube assembly chrome with rear mount hose barb, for N-EVAPs made after 1997</p> <p>Sold individually</p>	<p>P0627</p>	
<p>Swivel tube fitting for stainless steel gas supply tube</p> <p>Sold individually</p>	<p>P1204</p>	
<p>Paper filled air filter for N2 supply, Z-purge, or air compressor</p> <p>Sold individually</p>	<p>NA0403</p>	

Organomation®

Item Description	Part Number	Picture
Connector tube for gas supply, 6 ft coiled polyurethane plastic Sold individually	NA1101	
Gas supply tube Sold individually	P1230	
Hoist pulley and band spring Sold individually Please contact us before ordering to confirm correct band spring size	P1102 P1202 P1224 P1524	
Thumb Screw Sold individually	P0612-N	
Gas flow meter Sold individually	For 0-10 LPM: NA1221 For 0-20 LPM: NA1521 For 0-30 LPM: NA1421	
Pressure reducing regulator for gas, 30-150 psig in, 0-30 psig out	NA0630	

Organomation®

Item Description	Part Number	Picture
Thermometer, 6 in long, 100°C for N-EVAP water bath	NA1110	
Thermometer, 6 in long, 150°C for N-EVAP dry bath	NA1121	
Aluminum beads for dry bath	3.5 liters NA1119	
Sold in 3.5 liters for N-EVAP 111 and 5.5 liters for N-EVAP 112	5.5 liters NA1219	
Glass beads for acid resistant dry bath	3.5 liters NA1120	
Sold in 3.5 liters for N-EVAP 111 and 5.5 liters for N-EVAP 112	5.5 liters NA1220	
12 position sample holder for small glassware, 5-16 mm in outside diameter	NA1132	
Compatible with N-EVAP 111		
6 position sample holder for large glassware, 30-56 mm in outside diameter	NA1113	
Compatible with N-EVAP 111		
8 position sample holder for large glassware, 30-70 mm in outside diameter	NA1228	
Compatible with N-EVAP 112		