










# Organomation®

## ROT-X-TRACT Parts and Accessories List




*\*\*Please contact us for part numbers for acid resistant units\*\**

Item Description	Part Number	Picture
<p>Corning Accelerated One-Step extractor, 1000 mL, per position</p> <p>Compatible with the 8 position ROT-X-TRACT-LC</p>	GS3380	 A clear plastic extractor with a bulb and a stopcock assembly.
<p>J-shaped Snyder column</p> <p>Sold individually</p>	GP3441	 A clear plastic J-shaped column with a stopcock.
<p>Stopcock assembly; includes stopcock, red caps, and tubing</p> <p>Sold individually</p>	GA3447	 A stopcock assembly with a red cap and a blue handle.
<p>10 mL Concentrator Tube</p> <p>Sold individually</p>	GP2243	 A clear plastic concentrator tube with a stopcock and a red cap.

# Organomation®

Item Description	Part Number	Picture
<p>Y-shaped drain cup</p> <p>Sold individually</p>	<p>GP3442</p>	
<p>Membrane support screen</p> <p>Sold individually</p>	<p>GA3449</p>	
<p>Teflon sleeve for 24/40 ground glass joints</p> <p>Sold individually</p>	<p>GA2286-H</p>	
<p>Glassware set includes a 37 mm Soxhlet extractor, 250 mL round flask, and an Allihn condenser with 45/40 ground glass joint; per position</p> <p>Compatible with the 8 and 10 position ROT-X-TRACT-S</p>	<p>GS3137</p>	
<p>Pressure reducing regulator for water supply</p> <p>Compatible with all ROT-X-TRACT models</p>	<p>XA0631</p>	

# Organomation®

Item Description	Part Number	Picture
<p>Corning Membranes</p> <p>Sold in packs of 20 or 200</p>	<p>Pack of 20 GA3446-20</p> <p>Pack of 200 GA3446-200</p>	
<p>Boiling grids</p> <p>Sold in sets of 100</p>	<p>GA2245</p>	
<p>Recirculating Chiller; Supports up to 12 condensers</p> <p>Compatible with the 8 position ROT-X-TRACT-S and ROT-X-TRACT-LC</p>	<p>17760-C2T</p>	
<p>Recirculating Chiller; Supports up to 16 condensers</p> <p>Compatible with the 10 position ROT-X-TRACT-S</p>	<p>17760-C3T</p>	