











S-EVAP-KD & S-EVAP-RB Parts and Accessories List






Please contact us for part numbers for acid resistant units

Item Description	Part Number	Picture
<p>Buna O-ring set for water manifold; Includes one inner and one outer O-ring</p> <p>Please contact us for installation instructions</p>	P1332	
<p>White silicone tubing for water supply, 1/4 in ID x 3/8 in OD</p> <p>Sold in 10 ft or 25 ft sections</p>	<p>10 Feet: P2616</p> <p>25 Feet: P2613</p>	
<p>Precut white silicone tubing for water supply with elbow fitting</p> <p>Sold individually</p>	P1801	
<p>Bath covers for unused sample positions</p>	<p>2.125 inch diameter XA3283</p> <p>3.0 inch diameter XA2282</p> <p>3.875 inch diameter XA2283</p>	
<p>Thermometer, 6 in long, 100°C for analog baths</p>	NA1110	

Organomation®

Item Description	Part Number	Picture
<p>Glassware set for individual collection, per position</p> <p>Available for 250 mL or 500 mL KD flasks</p>	<p>250 mL KD flask set: GS2152</p> <p>500 mL KD flask set: GS2153</p>	
<p>Glassware set for central collection, per position</p> <p>Available for 250 mL or 500 mL KD flasks</p>	<p>250 mL KD flask set: GS2152-C</p> <p>500 mL KD flask set: GS2153-C</p>	
<p>Central collection kit</p> <p>Available for: 8 Position, 2.5 L for 12027 8 Position, 4 L for 12018 10 Position, 4 L for 12010</p>	<p>8 Position, 2.5 L GS2156-C8</p> <p>8 Position, 4 L GS2158-C8</p> <p>10 Position, 4 L GS2158-C10</p>	
<p>Rack to hold KD flasks when not in use. Compatible with 250 mL and 500 mL KD flasks.</p> <p>Sold individually</p>	<p>XA2290</p>	
<p>Nitrogen purge manifold to dry evaporation flasks</p> <p>Available for: 8 Position S-EVAP-RB 10 Position S-EVAP-RB</p>	<p>8 Position -N8</p> <p>10 Position -N10</p>	

Organomation®

Item Description	Part Number	Picture
<p>Glassware set for solvent collection, per position</p> <p>Available for 250 mL or 500 mL RB flasks</p>	<p>250 mL RB flask set: GS2162</p> <p>500 mL RB flask set: GS2163</p>	
<p>Adapter for gas injection, per position</p> <p>Necessary if choosing the nitrogen purge option</p>	GP3224	
<p>Stainless steel adapter ring for use with smaller glassware, per position</p> <p>Available for 125 mL flasks or 250 mL flasks</p>	<p>For 125 mL flasks: XA2286</p> <p>For 250 mL flasks: XA2284</p>	
<p>Pressure reducing regulator for water supply</p> <p>Sold individually</p>	XA0631	
<p>Boiling grids</p> <p>Sold in sets of 100</p>	GA2245	
<p>Teflon sleeve for 24/40 ground glass joints</p> <p>Sold individually</p>	GA2286-H	