











Organomation®

S-EVAP-KD & S-EVAP-RB Parts and Accessories List






Please contact us for part numbers for acid resistant units

Item Description	Part Number	Picture
Buna O-ring set for water manifold; Includes one inner and one outer O-ring Please contact us for installation instructions	P1332	
White silicone tubing for water supply, 1/4 in ID x 3/8 in OD Sold in 10 ft or 25 ft sections	10 Feet: P2616 25 Feet: P2613	
Precut white silicone tubing for water supply with elbow fitting Sold individually	P1801	
Bath covers for unused sample positions	2.125 inch diameter XA3283 3.0 inch diameter XA2282 3.875 inch diameter XA2283	
Thermometer, 6 in long, 100°C for analog baths	NA1110	

Organomation®

Item Description	Part Number	Picture
<p>Glassware set for individual collection, per position</p> <p>Available for 250 mL or 500 mL KD flasks</p>	<p>250 mL KD flask set: GS2152</p> <p>500 mL KD flask set: GS2153</p>	
<p>Glassware set for central collection, per position</p> <p>Available for 250 mL or 500 mL KD flasks</p>	<p>250 mL KD flask set: GS2152-C</p> <p>500 mL KD flask set: GS2153-C</p>	
<p>Central collection kit</p> <p>Available for: 8 Position, 2.5 L for 12027 8 Position, 4 L for 12018 10 Position, 4 L for 12010</p>	<p>8 Position, 2.5 L GS2156-C8</p> <p>8 Position, 4 L GS2158-C8</p> <p>10 Position, 4 L GS2158-C10</p>	
<p>Rack to hold KD flasks when not in use. Compatible with 250 mL and 500 mL KD flasks.</p> <p>Sold individually</p>	<p>XA2290</p>	
<p>Nitrogen purge manifold to dry evaporation flasks</p> <p>Available for: 8 Position S-EVAP-RB 10 Position S-EVAP-RB</p>	<p>8 Position -N8</p> <p>10 Position -N10</p>	

Organomation®

Item Description	Part Number	Picture
<p>Glassware set for solvent collection, per position</p> <p>Available for 250 mL or 500 mL RB flasks</p>	<p>250 mL RB flask set: GS2162</p> <p>500 mL RB flask set: GS2163</p>	
<p>Adapter for gas injection, per position</p> <p>Necessary if choosing the nitrogen purge option</p>	<p>GP3224</p>	
<p>Stainless steel adapter ring for use with smaller glassware, per position</p> <p>Available for 125 mL flasks or 250 mL flasks</p>	<p>For 125 mL flasks: XA2284</p> <p>For 250 mL flasks: XA2286</p>	
<p>Pressure reducing regulator for water supply</p> <p>Sold individually</p>	<p>XA0631</p>	
<p>Boiling grids</p> <p>Sold in sets of 100</p>	<p>GA2245</p>	
<p>Teflon sleeve for 24/40 ground glass joints</p> <p>Sold individually</p>	<p>GA2286-H</p>	