



APRIA®

Guidelines for Replacement Sleep Supplies

| ACCESSORIES | RE-ORDER |
|---------------------|----------------|
| Mask | 1 per 3 months |
| Nasal Pillow | 2 per 1 month |
| Headgear | 1 per 6 months |
| Chin Strap | 1 per 6 months |
| Tubing | 1 per 3 months |
| Disposable Filters | 2 per 1 month |
| Humidifier Chambers | 1 per 6 months |

**TO REORDER CPAP
SUPPLIES PLEASE CALL**

1-888-452-4363